

REEP5

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human REEP5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DP1; TB2; YOP1; D5S346; C5orf18;

To place an order, please [Click HERE](#).

Catalog #: A4191

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7905

Isotype: IgG

Swiss Prot: Q00765

Purity: Affinity purification

For research use only.