

STAM

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STAM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

STAM1; STAM-1;

Catalog #: A4198

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8027

Isotype: IgG

Swiss Prot: Q92783

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the signal-transducing adaptor molecule family. These proteins mediate downstream signaling of cytokine receptors and also play a role in ER to Golgi trafficking by interacting with the coat protein II complex. The encoded protein also associates with hepatocyte growth factor-regulated substrate to form the endosomal sorting complex required for transport-0 (ESCRT-0), which sorts ubiquitinated membrane proteins to the ESCRT-1 complex for lysosomal degradation. Alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).