

IFT88

Reactivity: Human Mouse Rat Dog

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:50

Calculated MW: 94kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human IFT88

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

d

Synonym:

DAF19; TG737; TTC10; hTg737; D13S1056E;

Catalog #: A4204

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8100

Isotype: IgG

Swiss Prot: Q13099

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the tetratrico peptide repeat (TPR) family. Mutations of a similar gene in mouse can cause polycystic kidney disease. Two transcript variants encoding distinct isoforms have been identified for this gene.

To place an order, please [Click HERE](#).