

USP11

Reactivity: Human

Tested applications: WB FC

Recommended Dilution: WB 1:500 - 1:2000 FC 1:20 - 1:50

Calculated MW: 110kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human USP11

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UHX1;

Catalog #: A4208

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 8237

Isotype: IgG

Swiss Prot: P51784

Purity: Affinity purification

For research use only.

Background:

Protein ubiquitination controls many intracellular processes, including cell cycle progression, transcriptional activation, and signal transduction. This dynamic process, involving ubiquitin conjugating enzymes and deubiquitinating enzymes, adds and removes ubiquitin. Deubiquitinating enzymes are cysteine proteases that specifically cleave ubiquitin from ubiquitin-conjugated protein substrates. This gene encodes a deubiquitinating enzyme which lies in a gene cluster on chromosome Xp11.

To place an order, please [Click HERE](#).