

RBM10

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 104kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RBM10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

S1-1; TARPS; GPATC9; ZRANB5; GPATCH9; DXS8237E;

Catalog #: A4209

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8241

Isotype: IgG

Swiss Prot: P98175

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein that belongs to a family proteins that contain an RNA-binding motif. The encoded protein associates with hnRNP proteins and may be involved in regulating alternative splicing. Defects in this gene are the cause of the X-linked recessive disorder, TARP syndrome. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).