

## FZD7

**Reactivity:** Human Mouse

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 64kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human FZD7

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

t

**Synonym:**

FzE3;

**Catalog #:** A4213

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 8324

**Isotype:** IgG

**Swiss Prot:** O75084

**Purity:** Affinity purification

For research use only.

**Background:**

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD7 protein contains an N-terminal signal sequence, 10 cysteine residues typical of the cysteine-rich extracellular domain of Fz family members, 7 putative transmembrane domains, and an intracellular C-terminal tail with a PDZ domain-binding motif. FZD7 gene expression may downregulate APC function and enhance beta-catenin-mediated signals in poorly differentiated human esophageal carcinomas.

*To place an order, please [Click HERE](#).*