

RABEP1

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 99kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RABEP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

RAB5EP; RABPT5;

To place an order, please [Click HERE](#).

Catalog #: A4273

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9135

Isotype: IgG

Swiss Prot: Q15276

Purity: Affinity purification

For research use only.