

VAPA

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human VAPA

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VAP-A; VAP33; VAP-33; hVAP-33;

Catalog #: A4279

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9218

Isotype: IgG

Swiss Prot: Q9P0L0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type IV membrane protein. It is present in the plasma membrane and intracellular vesicles. It may also be associated with the cytoskeleton. This protein may function in vesicle trafficking, membrane fusion, protein complex assembly and cell motility.

Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.

To place an order, please [Click HERE](#).