

UBE2L6

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human UBE2L6

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1mg/ml

Synonym:

RIG-B; UBCH8;

Catalog #: A4282

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9246

Isotype: IgG

Swiss Prot: O14933

Purity: Affinity purification

For research use only.

Background:

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is highly similar in primary structure to the enzyme encoded by the UBE2L3 gene. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

To place an order, please [Click HERE](#).