

PLAA

Reactivity: Human Mouse

Tested applications: WB IP

Recommended Dilution: WB 1:200 - 1:2000 IP 1:20 - 1:100

Calculated MW: 87kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PLAA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DOA1; PLAP; PLA2P;

To place an order, please [Click HERE](#).

Catalog #: A4294

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9373

Isotype: IgG

Swiss Prot: Q9Y263

Purity: Affinity purification

For research use only.