

EIF4E2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EIF4E2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

4EHP; IF4e; 4E-LP; EIF4EL3;

To place an order, please [Click HERE](#).

Catalog #: A4305

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9470

Isotype: IgG

Swiss Prot: O60573

Purity: Affinity purification

For research use only.