

RPH3AL

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RPH3AL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NOC2;

Catalog #: A4311

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9501

Isotype: IgG

Swiss Prot: Q9UNE2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene plays a direct regulatory role in calcium-ion-dependent exocytosis in both endocrine and exocrine cells and plays a key role in insulin secretion by pancreatic cells. This gene is likely a tumor suppressor. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).