

FXR2

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 74kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human FXR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

FXR2P; FMR1L2;

Catalog #: A4313

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9513

Isotype: IgG

Swiss Prot: P51116

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a RNA binding protein containing two KH domains and one RCG box, which is similar to FMRP and FXR1. It associates with polyribosomes, predominantly with 60S large ribosomal subunits. This encoded protein may self-associate or interact with FMRP and FXR1. It may have a role in the development of fragile X mental retardation syndrome.

To place an order, please [Click HERE](#).