

PREPL

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 84kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PREPL

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A4322

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9581

Isotype: IgG

Swiss Prot: Q4J6C6

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.

To place an order, please [Click HERE](#).