

MTRF1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human MTRF1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RF1; MRF1; MTTRF1;

Catalog #: A4327

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9617

Isotype: IgG

Swiss Prot: O75570

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was determined by in silico methods to be a mitochondrial protein with similarity to the peptide chain release factors (RFs) discovered in bacteria and yeast. The peptide chain release factors direct the termination of translation in response to the peptide chain termination codons. Initially thought to have a role in the termination of mitochondria protein synthesis, a recent publication found no mitochondrial translation release functionality. Multiple alternatively spliced transcript variants have been suggested by mRNA and EST data; however, their full-length natures are not clear.

To place an order, please [Click HERE](#).