

CLINT1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

FC 1:20 - 1:50

Calculated MW: 68kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human CLINT1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ENTH; EPN4; EPNR; CLINT;

Catalog #: A4332

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9685

Isotype: IgG

Swiss Prot: Q14677

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with similarity to the epsin family of endocytic adapter proteins. The encoded protein interacts with clathrin, the adapter protein AP-1 and phosphoinositides. This protein may be involved in the formation of clathrin coated vesicles and trafficking between the trans-Golgi network and endosomes. Mutations in this gene are associated with a susceptibility to schizophrenia and psychotic disorders. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).