

TOM70

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 67kDa

Observed MW: Refer to figures

Immunogen:

Recombinant Protein of human TOM70

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A4349

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9868

Isotype: IgG

Swiss Prot: O94826

Purity: Affinity purification

For research use only.

Background:

This gene encodes an import receptor of the outer mitochondrial membrane that is part of the translocase of the outer membrane complex. This protein is involved in the import of mitochondrial precursor proteins.

To place an order, please [Click HERE](#).