

GNPDA1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GNPDA1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GPI; HLN; GNP1; GNPI; GNPDA;

Catalog #: A4358

Antibody Type:

Monoclonal Antibody

Species: 23

Gene ID: 10007

Isotype: IgG

Swiss Prot: P46926

Purity: Affinity purification

For research use only.

Background:

Glucosamine-6-phosphate deaminase (EC 3.5.99.6) is an allosteric enzyme that catalyzes the reversible conversion of D-glucosamine-6-phosphate into D-fructose-6-phosphate and ammonium (Arreola et al., 2003)

To place an order, please [Click HERE](#).