

MAPKAP1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human MAPKAP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MIP1; SIN1; JC310; SIN1b; SIN1g;

Catalog #: A4960

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79109

Isotype: IgG

Swiss Prot: Q9BPZ7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is highly similar to the yeast SIN1 protein, a stress-activated protein kinase. Alternatively spliced transcript variants encoding distinct isoforms have been described. Alternate polyadenylation sites as well as alternate 3' UTRs have been identified for transcripts of this gene.

To place an order, please [Click HERE](#).