

TNIP2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:50 - 1:100

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TNIP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KLIP; ABIN2; FLIP1;

Catalog #: A4962

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79155

Isotype: IgG

Swiss Prot: Q8NFZ5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which acts as an inhibitor of NFkappaB activation. The encoded protein is also involved in MAP/ERK signaling pathway in specific cell types. It may be involved in apoptosis of endothelial cells. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.

To place an order, please [Click HERE](#).