

RBM42

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RBM42

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Catalog #: A4964

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79171

Isotype: IgG

Swiss Prot: Q9BTD8

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).