

GTDC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GTDC1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

mat-Xa; Hmat-Xa;

To place an order, please [Click HERE](#).

Catalog #: A4980

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79712

Isotype: IgG

Swiss Prot: Q4AE62

Purity: Affinity purification

For research use only.