

NHEJ1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NHEJ1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Synonym:

XLF;

Catalog #: A4985

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79840

Isotype: IgG

Swiss Prot: Q9H9Q4

Purity: Affinity purification

For research use only.

Background:

Double-strand breaks in DNA result from genotoxic stresses and are among the most damaging of DNA lesions. This gene encodes a DNA repair factor essential for the nonhomologous end-joining pathway, which preferentially mediates repair of double-stranded breaks. Mutations in this gene cause different kinds of severe combined immunodeficiency disorders.

To place an order, please [Click HERE](#).