

LPCAT1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LPCAT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AYTL2; lpcat; PFAAP3;

Catalog #: A4987

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79888

Isotype: IgG

Swiss Prot: Q8NF37

Purity: Affinity purification

For research use only.

Background:

Lysophosphatidylcholine (LPC) acyltransferase (LPCAT; EC 2.3.1.23) catalyzes the conversion of LPC to phosphatidylcholine (PC) in the remodeling pathway of PC biosynthesis (Nakanishi et al., 2006

To place an order, please [Click HERE](#).