

RMI1

Reactivity: Human

Tested applications: WB IHC IP

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RMI1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Synonym:

BLAP75; FAAP75; C9orf76;

Catalog #: A4991

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80010

Isotype: IgG

Swiss Prot: Q9H9A7

Purity: Affinity purification

For research use only.

Background:

RMI1 is a component of protein complexes that limit DNA crossover formation via the dissolution of double Holliday junctions (Raynard et al., 2006)

To place an order, please [Click HERE](#).