

GRPEL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GRPEL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

GrpE; HMGE;

To place an order, please [Click HERE](#).

Catalog #: A4999

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80273

Isotype: IgG

Swiss Prot: Q9HAV7

Purity: Affinity purification

For research use only.