

ZBP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 46kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DAI; DLM1; DLM-1; C20orf183;

Catalog #: A5011

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81030

Isotype: IgG

Swiss Prot: Q9H171

Purity: Affinity purification

For research use only.

Background:

This gene encodes a Z-DNA binding protein. The encoded protein plays a role in the innate immune response by binding to foreign DNA and inducing type-I interferon production.

Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).