

DOCK8

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 239kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human DOCK8

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MRD2; ZIR8; HEL-205;

Catalog #: A5021

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81704

Isotype: IgG

Swiss Prot: Q8NF50

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the DOCK180 family of guanine nucleotide exchange factors. Guanine nucleotide exchange factors interact with Rho GTPases and are components of intracellular signaling networks. Mutations in this gene result in the autosomal recessive form of the hyper-IgE syndrome. Alternatively spliced transcript variants encoding different isoforms have been described.

To place an order, please [Click HERE](#).