

EIF2A

Reactivity: Human

Tested applications: WB IHC FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 FC 1:20 - 1:50

Calculated MW: 65kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human EIF2A

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CDA02; EIF-2A; MST089; MSTP004; MSTP089;

Catalog #: A5035

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83939

Isotype: IgG

Swiss Prot: Q9BY44

Purity: Affinity purification

For research use only.

Background:

EIF2A is a 65-kD protein that catalyzes the formation of puromycin-sensitive 80S preinitiation complexes (Zoll et al., 2002)

To place an order, please [Click HERE](#).