

MND1

Reactivity: Human Mouse

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human MND1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GAJ;

Catalog #: A5038

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84057

Isotype: IgG

Swiss Prot: Q9BWT6

Purity: Affinity purification

For research use only.

Background:

The product of the MND1 gene associates with HOP2 (MIM 608665) to form a stable heterodimeric complex that binds DNA and stimulates the recombinase activity of RAD51 (MIM 179617) and DMC1 (MIM 602721) (Chi et al., 2007)

To place an order, please [Click HERE](#).