

ANTXR1

Reactivity: Human

Tested applications: WB IHC

Calculated MW: 564

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human ANTXR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A5045

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84168

Isotype: IgG

Swiss Prot: Q9H6X2

Purity: Affinity purification

For research use only.