

SARNP

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 24kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human SARNP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HCC1; THO1; CIP29; HSPC316;

Catalog #: A5054

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84324

Isotype: IgG

Swiss Prot: P82979

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is upregulated in response to various cytokines. The encoded protein may play a role in cell cycle progression. A translocation between this gene and the myeloid/lymphoid leukemia gene, resulting in expression of a chimeric protein, has been associated with acute myelomonocytic leukemia. Pseudogenes exist on chromosomes 7 and 8. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).