

PARD6B

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PAR6B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

PAR6B;

Catalog #: A5060

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84612

Isotype: IgG

Swiss Prot: Q9BYG5

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the PAR6 family and encodes a protein with a PSD95/Discs-large/ZO1 (PDZ) domain, an OPR domain and a semi-Cdc42/Rac interactive binding (CRIB) domain. This cytoplasmic protein is involved in asymmetrical cell division and cell polarization processes as a member of a multi-protein complex.

To place an order, please [Click HERE](#).