

ACCS

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human ACCS

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ACS; PHACS;

To place an order, please [Click HERE](#).

Catalog #: A5062

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84680

Isotype: IgG

Swiss Prot: Q96QU6

Purity: Affinity purification

For research use only.