

PIGO

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100

Calculated MW: 119kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PIGO

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HPMRS2;

Catalog #: A5065

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84720

Isotype: IgG

Swiss Prot: Q8TEQ8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is involved in glycosylphosphatidylinositol (GPI)-anchor biosynthesis. The GPI-anchor is a glycolipid which contains three mannose molecules in its core backbone. The GPI-anchor is found on many blood cells and serves to anchor proteins to the cell surface. This protein is involved in the transfer of ethanolaminephosphate (EtNP) to the third mannose in GPI. At least three alternatively spliced transcripts encoding two distinct isoforms have been found for this gene.

To place an order, please [Click HERE](#).