

GADD45GIP1

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:200 - 1:2000 IP 1:20 - 1:50

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GADD45GIP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRG6; CRIF1; PLINP; CKBBP2; Plinp1; MRP-L59; PLINP-1; CKbetaBP2;

Catalog #: A5081

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 90480

Isotype: IgG

Swiss Prot: Q8TAE8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear-localized protein that may be induced by p53 and regulates the cell cycle by inhibiting G1 to S phase progression. The encoded protein may interact with other cell cycle regulators.

To place an order, please [Click HERE](#).