

CREB3L1

Reactivity: Human Mouse

Tested applications: WB IHC FC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 FC:20 - 1:50

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human CREB3L1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OASIS;

Catalog #: A5082

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 90993

Isotype: IgG

Swiss Prot: Q96BA8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is normally found in the membrane of the endoplasmic reticulum (ER). However, upon stress to the ER, the encoded protein is cleaved and the released cytoplasmic transcription factor domain translocates to the nucleus. There it activates the transcription of target genes by binding to box-B elements.

To place an order, please [Click HERE](#).