

MOB1B

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human MOB1B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MATS2; MOB4A; MOBKL1A;

Catalog #: A5091

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92597

Isotype: IgG

Swiss Prot: Q7L9L4

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).