

HIST3H2A

Reactivity: Human Mouse

Tested applications: WB

Calculated MW: 130

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human HIST3H2A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A5093

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92815

Isotype: IgG

Swiss Prot: Q7L7L0

Purity: Affinity purification

For research use only.