

IMP4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human IMP4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BXDC4;

Catalog #: A5095

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92856

Isotype: IgG

Swiss Prot: Q96G21

Purity: Affinity purification

For research use only.

Background:

IMP4 forms a ternary complex with IMP3 (MIM 612980) and MPP10 (MPHOSPHO10; MIM 605503) that interacts with U3 small nucleolar RNA (snoRNA), which is required for the early cleavage steps in pre-rRNA processing (Granneman et al., 2003)

To place an order, please [Click HERE](#).