

PRKCDBP

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PRKCDBP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SRBC; HSRBC; CAVIN3; cavin-3;

Catalog #: A5105

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 112464

Isotype: IgG

Swiss Prot: Q969G5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was identified as a binding protein of the protein kinase C, delta (PRKCD). The expression of this gene in cultured cell lines is strongly induced by serum starvation. The expression of this protein was found to be down-regulated in various cancer cell lines, suggesting the possible tumor suppressor function of this protein.

To place an order, please [Click HERE](#).