

PIK3IP1

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PIK3IP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HGFL; hHGFL(S);

To place an order, please [Click HERE](#).

Catalog #: A5110

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 113791

Isotype: IgG

Swiss Prot: Q96FE7

Purity: Affinity purification

For research use only.