

KCTD12

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:10 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human KCTD12

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

g

Synonym:

PFET1; PFETIN; C13orf2;

To place an order, please [Click HERE](#).

Catalog #: A5116

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 115207

Isotype: IgG

Swiss Prot: Q96CX2

Purity: Affinity purification

For research use only.