

ACAP3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 92kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human ACAP3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CENTB5;

To place an order, please [Click HERE](#).

Catalog #: A5123

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116983

Isotype: IgG

Swiss Prot: Q96P50

Purity: Affinity purification

For research use only.