

NEDD1

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NEDD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

GCP-WD; TUBGCP7;

To place an order, please [Click HERE](#).

Catalog #: A5128

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 121441

Isotype: IgG

Swiss Prot: Q8NHV4

Purity: Affinity purification

For research use only.