

NEU4

Reactivity: Human

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NEU4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5141

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 129807

Isotype: IgG

Swiss Prot: Q8WWR8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to a family of glycohydrolytic enzymes, which remove terminal sialic acid residues from various sialo derivatives, such as glycoproteins, glycolipids, oligosaccharides, and gangliosides. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

To place an order, please [Click HERE](#).