

PIP5KL1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PIP5KL1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PIPKH; bA203J24.5;

Catalog #: A5152

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 138429

Isotype: IgG

Swiss Prot: Q5T9C9

Purity: Affinity purification

For research use only.

Background:

PIP5KL1 is a phosphoinositide kinase-like protein that lacks intrinsic lipid kinase activity but associates with type I PIPKs (see PIP5K1A; MIM 603275) and may play a role in localization of PIPK activity (Chang et al., 2004)

To place an order, please [Click HERE](#).