

ASB9

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human ASB9

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5156

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 140462

Isotype: IgG

Swiss Prot: Q96DX5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ankyrin repeat and suppressor of cytokine signaling (SOCS) box protein family. Members of this family can interact with the elongin B-C adapter complex via their SOCS box domain and further complex with the cullin and ring box proteins to form E3 ubiquitin ligase complexes. They may function to mediate the substrate-recognition of the E3 ubiquitin ligases. A transcribed pseudogene of this gene has been identified on chromosome 15. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).