

TRIM69

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TRIM69

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bl

Synonym:

Trif; HSD34; RNF36; HSD-34;

Catalog #: A5157

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 140691

Isotype: IgG

Swiss Prot: Q86WT6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the RING-B-box-coiled-coil (RBCC) family and encodes a protein with an N-terminal RING finger motif, a PRY domain and a C-terminal SPRY domain. The mouse ortholog of this gene is specifically expressed in germ cells at the round spermatid stages during spermatogenesis and, when overexpressed, induces apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been described.

To place an order, please [Click HERE](#).