

VTI1A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VTI1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MMDS3; MVti1; VTI1RP2; Vti1-rp2;

To place an order, please [Click HERE](#).

Catalog #: A5162

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 143187

Isotype: IgG

Swiss Prot: Q96AJ9

Purity: Affinity purification

For research use only.